

QUADRASTEP SYSTEM®

THANK YOU For purchasing our
QUADRASTEP® Starter Kit



Now let's help you get started!
Here is how you set up and use your Fit Kit:

HOW IS THE KIT SET UP?

Your **QUADRASTEP SYSTEM®** Fit Kit has 6 drawers, one for each size orthotic, and one for extra materials. You will see a label on each drawer indicating the size of the orthotics that should be stored in that drawer.

In the bottom drawer you will find your marketing materials in a portfolio, along with your waiting room brochure holder filled with **QUADRASTEP®** patient's brochures. These materials are enough to get you started, additional copies can be downloaded to print on our website nolaro24.com/kit.html. Look for the "marketing material" link.

WHAT WILL I FIND IN EACH DRAWER?

Your **QUADRASTEP SYSTEM®** Fit Kit has 6 drawers, one for each size **QUADRASTEP®** orthotic.

In the first drawer you will find Size 0 of the **E+ Orthotic (Dark Red)**, and size 1 of all of the orthotics, one each **QUAD A (Yellow)**, **QUAD B (Purple)**, **QUAD C (Blue)**, **QUAD D (Green)**, **QUAD E (Red)**, **E+ (Dark Red)** and **QUAD F (Orange)**. Drawer 2 will have all of the Quads (colors) in Size 2, and so on through the 5 drawers/5 sizes. In the 6th drawer you will find 2 pairs of orthotics, Quads **C** and **D** in size 6, as well as your marketing materials and brochure holder.

It is best to keep your orthotics organized in this way for easy access when fitting a patient. That way you will always know what is available and have it handy. Another trick is to always arrange your drawers so that the most commonly dispensed orthotics are in the front, and the least dispensed in the back, as shown.



SHOULD I SELL OUT OF MY KIT?

Your **QUADRASTEP SYSTEM®** Kit has one pair of orthotics in each size and each Quad, for a total of 38 pairs. We recommend that you keep your Kit available for fitting, then order a new pair of **QUADRASTEP®** orthotics to dispense to your patient. The orthotics in your Kit can be disinfected with any household disinfectant. You can also purchase a **QUADRASTEP SYSTEM®** Fit Kit for dispensing.